



# CALLISTO LARYNGOSCOPE RANGE

READY TO USE • FIBRE OPTIC • SINGLE USE • EXCEPTIONAL ILLUMINATION • LED • SETS

**TIMESCO**  
INSTRUMENTS FOR LIFE

# ABOUT THE RANGE

Callisto is the U.K's No.1 single use laryngoscope system. It has been designed to offer a true low cost alternative to reusable blades and handles with exceptional specifications to reduce cross infection, cost, maintenance and processing costs.

## PRODUCT FEATURES

### METAL SPATULA

Rigid metal spatula looks and feels like a reusable fibre optic blade.

### STRONG & SECURE

Metal hook bar offering exceptional strength for a single use handle and contains 2 metal ball bearings for secure fitting.

### EXCEPTIONAL ILLUMINATION

Highly polished 4mm acrylic core for maximum illumination.

### HIGH PERFORMANCE

Powerful bright LED bulb for excellent illumination with an average minimum output of 11,000 lux when measured 2cm from the bulb.

### EXTENSIVE RANGE

Available in a comprehensive range of convenient ready to use sets, including Eclipse (tilting tip) for elevation of the epiglottis; to provide an improved view of the larynx without excessive force.



### READY TO USE

Preloaded with alkaline button cell batteries, which are quick and easy to remove for safe battery disposal.



### DISPOSABLE

Features a low volume plastic construction to reduce plastic waste and cost of disposal.

## PRODUCT SPECIFICATION



CE marked



ISO 7376 compatible



Fibre-optic, green fitting



Single Use



Supplied with multilingual instructions



Latex Free

## RANGE INCLUDES



Blades



Handles



Sets

Supplied clinically clean and ready to use with a shelf life of:  
Blades: 5 years  
Handles: 2 years battery life

# LARYNGOSCOPE BLADES




## MACINTOSH



| Product Code   | Description                 | Measurements* | Quantity  |
|----------------|-----------------------------|---------------|-----------|
| DS.3940.150.05 | Callisto Macintosh Size 0   | 80mm x 13mm   | Box of 10 |
| DS.3940.150.10 | Callisto Macintosh Size 1   | 90mm x 18mm   | Box of 10 |
| DS.3940.150.15 | Callisto Macintosh Size 2   | 100mm x 18mm  | Box of 10 |
| DS.3940.150.20 | Callisto Macintosh Size 3   | 130mm x 18mm  | Box of 10 |
| DS.3940.150.23 | Callisto Macintosh Size 3.5 | 144mm x 18mm  | Box of 10 |
| DS.3940.150.25 | Callisto Macintosh Size 4   | 153mm x 18mm  | Box of 10 |
| DS.3940.150.30 | Callisto Macintosh Size 5   | 174mm x 16mm  | Box of 10 |

Macintosh design features a curved blade, are inserted into the vallecula, epiglottis fold, and moved laterally to expose the laryngeal cavity and vocal cords. These are the most popular blades worldwide.


## MILLER



| Product Code   | Description             | Measurements* | Quantity  |
|----------------|-------------------------|---------------|-----------|
| DS.3940.185.03 | Callisto Miller Size 00 | 65mm x 12mm   | Box of 10 |
| DS.3940.185.05 | Callisto Miller Size 0  | 79mm x 12mm   | Box of 10 |
| DS.3940.185.10 | Callisto Miller Size 1  | 105mm x 13mm  | Box of 10 |
| DS.3940.185.15 | Callisto Miller Size 2  | 158mm x 15mm  | Box of 10 |
| DS.3940.185.20 | Callisto Miller Size 3  | 197mm x 16mm  | Box of 10 |
| DS.3940.185.25 | Callisto Miller Size 4  | 208mm x 18mm  | Box of 10 |

Miller designs feature a straight blade, are inserted over the epiglottis and moved laterally to expose the laryngeal cavity and vocal cords. The smaller sizes 00, 0 and 1 are preferred by clinicians for paediatric and infant use.


## ROBERTSHAW



| Product Code   | Description                | Measurements* | Quantity  |
|----------------|----------------------------|---------------|-----------|
| DS.3940.235.05 | Callisto Robertshaw Size 0 | 94mm x 10mm   | Box of 10 |
| DS.3940.235.10 | Callisto Robertshaw Size 1 | 112mm x 12mm  | Box of 10 |

Due to the low profile on the Robertshaw blades, these are preferred by clinicians for paediatric and infant use.

## SEWARD



| Product Code   | Description            | Measurements* | Quantity  |
|----------------|------------------------|---------------|-----------|
| DS.3940.240.10 | Callisto Seward Size 1 | 107mm x 13mm  | Box of 10 |

Due to the low profile on the Seward blades, these are preferred by clinicians for paediatric and infant use.


In this brochure, only our Callisto single use fibre optic light guide laryngoscope systems are illustrated. We also manufacture the reusable laryngoscope systems, [Optima CLX](#), with removable fibre optic light cores, [Sirius](#) with fixed fibre optic light cores and Orion our standard light systems. Please contact us or refer to our website & separate publications for more information on these systems.

# LARYNGOSCOPE BLADES

Highly polished 4mm acrylic core for maximum illumination.

Supplied in boxes of 10 when sold individually or boxes of 5 when sold as a set.

Rigid metal spatula looks and feels like a reusable fibre optic blade.



No touch hook prevents the blade from touching the handle.

Green ISO 7376 fitting plastic block compatible with green fitting single use and reusable handles.

Single use logo clearly marked on the product.

2 metal ball bearings for secure fitting.



Tilting tip for elevation of the epiglottis, providing an improved view of the larynx, without excessive force.

## CALLISTO ECLIPSE BLADES

| Product Code   | Description                   | Measurements* | Quantity  |
|----------------|-------------------------------|---------------|-----------|
| DS.3940.175.20 | Callisto Eclipse Blade Size 3 | 130mm x 18mm  | Box of 10 |
| DS.3940.175.25 | Callisto Eclipse Blade Size 4 | 150mm x 18mm  | Box of 10 |

# LARYNGOSCOPE HANDLES

Metal hook bar offering exceptional strength for a single use handle.

This handle is ISO 7376 fitting and compatible with green fitting blades.

Preloaded with alkaline button cell batteries, which are quick and easy to remove for safe battery disposal.

Powerful, bright LED bulb for excellent illumination with a average minimum output of 11,000 lux when measured 2cm from the bulb.

Supplied in boxes of 25.

Low volume plastic construction reducing plastic waste and cost of disposal.

Clearly marked as single use on the base of the handle.



Medium



Compact

| Product Code   | Description                        | Measurements* | Quantity  |
|----------------|------------------------------------|---------------|-----------|
| DS.3940.480.10 | Callisto Preloaded Handle, Medium  | 139mm x 33mm  | Box of 25 |
| DS.3940.460.10 | Callisto Preloaded Handle, Compact | 120mm x 29mm  | Box of 25 |

## CALLISTO BLADE & LED PRELOADED MEDIUM HANDLE SETS



### CALLISTO MACINTOSH BLADE & LED PRELOADED MEDIUM HANDLE SET

| Product Code   | Description                                                | Quantity  |
|----------------|------------------------------------------------------------|-----------|
| DS.3950.950.05 | Callisto Blade Macintosh 0 & LED Preloaded Medium Handle   | Box of 10 |
| DS.3950.950.10 | Callisto Blade Macintosh 1 & LED Preloaded Medium Handle   | Box of 10 |
| DS.3950.950.15 | Callisto Blade Macintosh 2 & LED Preloaded Medium Handle   | Box of 10 |
| DS.3950.950.20 | Callisto Blade Macintosh 3 & LED Preloaded Medium Handle   | Box of 10 |
| DS.3950.950.23 | Callisto Blade Macintosh 3.5 & LED Preloaded Medium Handle | Box of 10 |
| DS.3950.950.25 | Callisto Blade Macintosh 4 & LED Preloaded Medium Handle   | Box of 10 |
| DS.3950.950.27 | Callisto Blade Macintosh 5 & LED Preloaded Medium Handle   | Box of 10 |



### CALLISTO MILLER BLADE & LED PRELOADED MEDIUM HANDLE SET

| Product Code   | Description                                            | Quantity  |
|----------------|--------------------------------------------------------|-----------|
| DS.3950.950.28 | Callisto Blade Miller 00 & LED Preloaded Medium Handle | Box of 10 |
| DS.3950.950.30 | Callisto Blade Miller 0 & LED Preloaded Medium Handle  | Box of 10 |
| DS.3950.950.35 | Callisto Blade Miller 1 & LED Preloaded Medium Handle  | Box of 10 |
| DS.3950.950.40 | Callisto Blade Miller 2 & LED Preloaded Medium Handle  | Box of 10 |
| DS.3950.950.45 | Callisto Blade Miller 3 & LED Preloaded Medium Handle  | Box of 10 |



### CALLISTO ECLIPSE (TILTING TIP) BLADE & LED PRELOADED MEDIUM HANDLE SET

| Product Code   | Description                                                      | Quantity |
|----------------|------------------------------------------------------------------|----------|
| DS.3950.953.20 | Callisto Blade Eclipse Macintosh 3 & LED Preloaded Medium Handle | Box of 5 |
| DS.3950.953.25 | Callisto Blade Eclipse Macintosh 4 & LED Preloaded Medium Handle | Box of 5 |



## CALLISTO BLADE & LED PRELOADED COMPACT HANDLE SETS



### CALLISTO MACINTOSH BLADE & LED PRELOADED COMPACT HANDLE SET

| Product Code   | Description                                                 | Quantity  |
|----------------|-------------------------------------------------------------|-----------|
| DS.3950.952.05 | Callisto Blade Macintosh 0 & LED Preloaded Compact Handle   | Box of 10 |
| DS.3950.952.10 | Callisto Blade Macintosh 1 & LED Preloaded Compact Handle   | Box of 10 |
| DS.3950.952.15 | Callisto Blade Macintosh 2 & LED Preloaded Compact Handle   | Box of 10 |
| DS.3950.952.20 | Callisto Blade Macintosh 3 & LED Preloaded Compact Handle   | Box of 10 |
| DS.3950.952.23 | Callisto Blade Macintosh 3.5 & LED Preloaded Compact Handle | Box of 10 |
| DS.3950.952.25 | Callisto Blade Macintosh 4 & LED Preloaded Compact Handle   | Box of 10 |
| DS.3950.952.27 | Callisto Blade Macintosh 5 & LED Preloaded Compact Handle   | Box of 10 |



### CALLISTO MILLER BLADE & LED PRELOADED COMPACT HANDLE SET

| Product Code   | Description                                             | Quantity  |
|----------------|---------------------------------------------------------|-----------|
| DS.3950.952.28 | Callisto Blade Miller 00 & LED Preloaded Compact Handle | Box of 10 |
| DS.3950.952.30 | Callisto Blade Miller 0 & LED Preloaded Compact Handle  | Box of 10 |
| DS.3950.952.35 | Callisto Blade Miller 1 & LED Preloaded Compact Handle  | Box of 10 |
| DS.3950.952.40 | Callisto Blade Miller 2 & LED Preloaded Compact Handle  | Box of 10 |
| DS.3950.952.45 | Callisto Blade Miller 3 & LED Preloaded Compact Handle  | Box of 10 |

## BENEFITS OF SINGLE USE

The Timesco Callisto single use sets have exceptional specifications to reduce cross infection at a cost low enough to justify single use. The sets contain single use blades and handles (preloaded with batteries) providing a ready to use disposable system to increase time efficiency and convenience.



Please refer to our website for our frequently asked questions and answers on the family of Callisto products. Visit: [www.timesco.com/frequently-asked-questions-single-use/](http://www.timesco.com/frequently-asked-questions-single-use/)



# PRODUCT LIST

| Product Code                                   | Description                        | Dimensions   | Quantity  |
|------------------------------------------------|------------------------------------|--------------|-----------|
| <b>CALLISTO MACINTOSH LARYNGOSCOPE BLADES</b>  |                                    |              |           |
| DS.3940.150.05                                 | Callisto Macintosh Size 0          | 80mm x 13mm  | Box of 10 |
| DS.3940.150.10                                 | Callisto Macintosh Size 1          | 90mm x 18mm  | Box of 10 |
| DS.3940.150.15                                 | Callisto Macintosh Size 2          | 100mm x 18mm | Box of 10 |
| DS.3940.150.20                                 | Callisto Macintosh Size 3          | 130mm x 18mm | Box of 10 |
| DS.3940.150.23                                 | Callisto Macintosh Size 3.5        | 144mm x 18mm | Box of 10 |
| DS.3940.150.25                                 | Callisto Macintosh Size 4          | 153mm x 18mm | Box of 10 |
| DS.3940.150.30                                 | Callisto Macintosh Size 5          | 174mm x 16mm | Box of 10 |
| <b>CALLISTO MILLER LARYNGOSCOPE BLADES</b>     |                                    |              |           |
| DS.3940.185.03                                 | Callisto Miller Size 00            | 65mm x 12mm  | Box of 10 |
| DS.3940.185.05                                 | Callisto Miller Size 0             | 79mm x 12mm  | Box of 10 |
| DS.3940.185.10                                 | Callisto Miller Size 1             | 105mm x 13mm | Box of 10 |
| DS.3940.185.15                                 | Callisto Miller Size 2             | 158mm x 15mm | Box of 10 |
| DS.3940.185.20                                 | Callisto Miller Size 3             | 197mm x 16mm | Box of 10 |
| DS.3940.185.25                                 | Callisto Miller Size 4             | 208mm x 18mm | Box of 10 |
| <b>CALLISTO ROBERTSHAW LARYNGOSCOPE BLADES</b> |                                    |              |           |
| DS.3940.235.05                                 | Callisto Robertshaw Size 0         | 94mm x 10mm  | Box of 10 |
| DS.3940.235.10                                 | Callisto Robertshaw Size 1         | 112mm x 12mm | Box of 10 |
| <b>CALLISTO SEWARD LARYNGOSCOPE BLADES</b>     |                                    |              |           |
| DS.3940.240.10                                 | Callisto Seward Size 1             | 107mm x 13mm | Box of 10 |
| <b>CALLISTO ECLIPSE LARYNGOSCOPE BLADES</b>    |                                    |              |           |
| DS.3940.175.20                                 | Callisto Eclipse Blade Size 3      | 130mm x 18mm | Box of 10 |
| DS.3940.175.25                                 | Callisto Eclipse Blade Size 4      | 150mm x 18mm | Box of 10 |
| <b>CALLISTO LARYNGOSCOPE HANDLES</b>           |                                    |              |           |
| DS.3940.480.10                                 | Callisto Preloaded Handle, Medium  | 139mm x 33mm | Box of 25 |
| DS.3940.460.10                                 | Callisto Preloaded Handle, Compact | 120mm x 29mm | Box of 25 |

# PRODUCT LIST (CONTINUED)

| Product Code                                                                      | Description                                                      | Quantity  |
|-----------------------------------------------------------------------------------|------------------------------------------------------------------|-----------|
| <b>CALLISTO MACINTOSH BLADE &amp; LED PRELOADED MEDIUM HANDLE SET</b>             |                                                                  |           |
| DS.3950.950.05                                                                    | Callisto Blade Macintosh 0 & LED Preloaded Medium Handle         | Box of 10 |
| DS.3950.950.10                                                                    | Callisto Blade Macintosh 1 & LED Preloaded Medium Handle         | Box of 10 |
| DS.3950.950.15                                                                    | Callisto Blade Macintosh 2 & LED Preloaded Medium Handle         | Box of 10 |
| DS.3950.950.20                                                                    | Callisto Blade Macintosh 3 & LED Preloaded Medium Handle         | Box of 10 |
| DS.3950.950.23                                                                    | Callisto Blade Macintosh 3.5 & LED Preloaded Medium Handle       | Box of 10 |
| DS.3950.950.25                                                                    | Callisto Blade Macintosh 4 & LED Preloaded Medium Handle         | Box of 10 |
| DS.3950.950.27                                                                    | Callisto Blade Macintosh 5 & LED Preloaded Medium Handle         | Box of 10 |
| <b>CALLISTO MILLER BLADE &amp; LED PRELOADED MEDIUM HANDLE SET</b>                |                                                                  |           |
| DS.3950.950.28                                                                    | Callisto Blade Miller 00 & LED Preloaded Medium Handle           | Box of 10 |
| DS.3950.950.30                                                                    | Callisto Blade Miller 0 & LED Preloaded Medium Handle            | Box of 10 |
| DS.3950.950.35                                                                    | Callisto Blade Miller 1 & LED Preloaded Medium Handle            | Box of 10 |
| DS.3950.950.40                                                                    | Callisto Blade Miller 2 & LED Preloaded Medium Handle            | Box of 10 |
| DS.3950.950.45                                                                    | Callisto Blade Miller 3 & LED Preloaded Medium Handle            | Box of 10 |
| <b>CALLISTO ECLIPSE (TILTING TIP) BLADE &amp; LED PRELOADED MEDIUM HANDLE SET</b> |                                                                  |           |
| DS.3950.953.20                                                                    | Callisto Blade Eclipse Macintosh 3 & LED Preloaded Medium Handle | Box of 5  |
| DS.3950.953.25                                                                    | Callisto Blade Eclipse Macintosh 4 & LED Preloaded Medium Handle | Box of 5  |
| <b>CALLISTO MACINTOSH BLADE &amp; LED PRELOADED COMPACT HANDLE SET</b>            |                                                                  |           |
| DS.3950.952.05                                                                    | Callisto Blade Macintosh 0 & LED Preloaded Compact Handle        | Box of 10 |
| DS.3950.952.10                                                                    | Callisto Blade Macintosh 1 & LED Preloaded Compact Handle        | Box of 10 |
| DS.3950.952.15                                                                    | Callisto Blade Macintosh 2 & LED Preloaded Compact Handle        | Box of 10 |
| DS.3950.952.20                                                                    | Callisto Blade Macintosh 3 & LED Preloaded Compact Handle        | Box of 10 |
| DS.3950.952.23                                                                    | Callisto Blade Macintosh 3.5 & LED Preloaded Compact Handle      | Box of 10 |
| DS.3950.952.25                                                                    | Callisto Blade Macintosh 4 & LED Preloaded Compact Handle        | Box of 10 |
| DS.3950.952.27                                                                    | Callisto Blade Macintosh 5 & LED Preloaded Compact Handle        | Box of 10 |
| <b>CALLISTO MILLER BLADE &amp; LED PRELOADED COMPACT HANDLE SET</b>               |                                                                  |           |
| DS.3950.952.28                                                                    | Callisto Blade Miller 00 & LED Preloaded Compact Handle          | Box of 10 |
| DS.3950.952.30                                                                    | Callisto Blade Miller 0 & LED Preloaded Compact Handle           | Box of 10 |
| DS.3950.952.35                                                                    | Callisto Blade Miller 1 & LED Preloaded Compact Handle           | Box of 10 |
| DS.3950.952.40                                                                    | Callisto Blade Miller 2 & LED Preloaded Compact Handle           | Box of 10 |
| DS.3950.952.45                                                                    | Callisto Blade Miller 3 & LED Preloaded Compact Handle           | Box of 10 |

# Contact Us

find out more today

For over 55 years' Timesco have been at the forefront of Medical product design and innovation. As a manufacturer and supplier of laryngoscope products, we have an unrivalled portfolio supplying both the UK directly and internationally through distributors.

As a family owned business, with headquarters in Basildon England, we strive to deliver products made to an exceptionally high standard; exceeding internationally recognised quality standards and customer expectations. We are recognised as one the most progressive and innovative companies within the industry and have a dedicated product development team which focuses on the changing needs of our clients and advances in clinical practice.

Contact Timesco today for further information or to request product samples.

**t** +44 (0)1268 297710    **e** [sales@timesco.com](mailto:sales@timesco.com)

**f** +44 (0)1268 297801    **w** [timesco.com](http://timesco.com)



ISO 13485



CE 1639



FDA Manufacturers  
License (USA)  
(Registration Number: 8021788)



SFDA  
CHINA



Ministry of  
Health

Classification: TIMESCO products are classified under the Harmonised Tariff: 90-1890-8400.

Timesco Healthcare Ltd Timesco House 3 Carnival Park  
Carnival Close Basildon Essex SS14 3WN England

**T** +44 (0) 1268 297 700 | **F** +44 (0) 1268 297 800

**E** [info@timesco.com](mailto:info@timesco.com) | [timesco.com](http://timesco.com)

**TIMESCO**  
INSTRUMENTS FOR LIFE